

# TIC Titanium Metal Powder Filter Cartridge

Cobetter **TIC**<sup>®</sup> Titanium Metal Powder Filter Cartridges composed of high-purity industrial-grade titanium powder (99.4%) with all elements sintered at high temperatures. Its features include anti-chemical corrosion, oxidation and high temperature resistance, and long service life.

As it is a low viscosity liquid filter, this filter results in good solid-liquid separation efficiency.

This filter is mainly used as a chemical filter to remove ozone-depleting substance and for the removal of carbon dioxide in food, pharmaceutical, and water treatment applications.

## Features and Benefits



- High-purity titanium construction
- Anti-corrosive; high temperature and oxidation resistant
- Uniform structure with narrow pore size distribution and high filtration efficiency
- No free-falling particles
- High porosity, low filtration resistance and high filtration efficiency
- Good compatibility with human tissue and blood due to its non-toxic and non-magnetic nature



## Materials of Construction(Five Layers)

Filter Layer	High-purity Titanium
End Cap	High-purity Titanium
Screw Cap	304 Stainless Steel
Reinforcing Layers	304/316 Stainless Steel

## Configurations

- Double Open-End (DOE)
- Single Open-End (SOE)

## Nominal Dimensions


Diameters	60mm
-----------	------

• Additional Diameter Specifications Available Upon Request

## Operating Conditions

- Max. Differential Pressure 3.0 bar / 21°C (forward flow)
- Max. Operating Temperature 280°C

## Manufacturing Process of TIC Titanium Metal Powder Filter Cartridges



Trap/Pre-Filtration

Microbiological Stabilization

Gas Filtration

Additional Filters



## Parameters

Code	Liquid Pore Size (µm)	Removal Ratings(µm)	Pore Efficiency	Absolute Removal Rating (µm) ②	Average Air Permeability (L/dm²min) ③	Flow Rate (m³/h) ④
1	0.45	32	30-50%	6	0.02	0.18
2	1.0	25		10	0.1	0.27
3	3.0	6.1		20	0.5	0.33
4	5.0	3.2		30	1.1	1.32
5	10	3.0		50	2.7	4.2
6	20	2.8		70	5.6	5.6
7	30	-		-	6.5	-
8	50	-		-	10.5	-
9	80	-		-	14.9	-
10	100	-		-	18	-
11	120	-		-	20	-

② Bubble Point Testing

③ Tested according to GB/T8786; Differential Pressure of 200Pa (in air)

④ Liquid Viscosity of 1 CP-S; diameter of 65mm; length of 10inches; pressure of 1.0bar

## Particle Efficiency

Particle Range	0.45µm	1 µm	3 µm	5 µm	10 µm
≥2µm	99.916%	99.895%	99.769%	82.546%	82.371%
≥5µm	99.974%	99.965%	99.910%	96.283%	96.079%
≥10µm	99.990%	99.986%	99.973%	98.875%	98.902%
≥12µm	99.987%	99.987%	99.986%	98.998%	98.982%
≥25µm	100.000%	100.000%	100.000%	99.996%	99.916%
≥35µm	100.000%	100.000%	100.000%	100.000%	99.966%
≥50µm	100.000%	100.000%	100.000%	100.000%	100.000%

## Ordering Information

TIC	Removal Ratings	End Cap	Nominal Length	Diameter	Seal Material	-F
0045=0.45µm	3000 =30µm	DOE =Double open end	05 = 5"	D22 =22mm	D70 =70mm	S =Silicon
0100=1.0µm	5000 =50µm	TC =222/Flat	10 =10"	D30 =30mm	D75 =75mm	E =EPDM
0300=3.0µm	8000 =80µm	SC =226/Flat	20 =20"	D40 =40mm	D80 =80mm	V =Viton
0500=5.0µm	100H =100µm	L =Screw	30 =30"	D50 =50mm	D120 =120mm	P =PFA/Viton
1000=10µm	120H =120µm		40 =40"	D60 =60mm		F =PTFE
2000=20µm						

## Cleaning Methods

Physical Cleaning Methods: Reverse-Flow by Clean Water; Reverse-Blow by Clear Air and Ultrasonic Wave

Chemical Cleaning Methods: Use Cleaning Agent Such As Diluted Acid, Diluted Alkalis, Oxidizer, and Surfactant

Contamination Materials	General Cleaning Method Procedures
Decarburization in Pharmaceutical and Chemical Industries	Reverse-blow and reverse-flow used more frequently; ultrasonic wave cleaning used when necessary
Non -Water Soluble Salts and Oxides in the Pharmaceutical Industry	Soak in 5% Concentration of Nitric Acid Solution
Original Liquid Filtration	Choose the correct cleaning methods as per the chemical properties of the contamination material; The Ultrasonic Wave Cleaning can be combined to use when necessary
Chemical Cleaning Methods	Detailed Procedures
Alkaline Cleaning	Alkaline Cleaning Soak filter in 3-5% Concentration of AR grade NaOH Solution for 30-60 minutes; solution temperature is 40°C. Flush the soaked filter inside, out with DI water or WFI water until the flushed solution turns neutral , and then test its conductivity. Dry with Pure Air ≥0.4Mpa
Acid Cleaning	Soak it in the 5% Concentration of Nitric Acid Solution for at least 8 hours; solution temperature is 40°C. Flush the soaked filter inside, out with DI water or WFI water until the flushed solution turns neutral, and then test its conductivity. Dry with Pure Air ≥ 0.4Mpa
Original Liquid Filtration	Clean filter with surfactant caused by contamination with Organic Pollution ( high concentration of Citric Acid recommended for Food and Beverage Applications)

## Length and Area<sup>®</sup>

Length	Filtration Area <sup>®</sup>
5 in. (125 mm)	0.024 m²
10 in. (300 mm)	0.056 m²
20 in. (500 mm)	0.094 m²
30 in. (750 mm)	0.141 m²
40 in. (1000 mm)	0.188 m²

⑤ Length and Other Sizes Are Customizable

⑥ Tested Filter Diameter is 65mm

