


H-CFH TFF Stainless Steel Holder



0.1 m² Stainless Steel Holder

Process development and small-volume parentals production with an EFA of 0.1m².

Unique design is fully compliant with biopharmaceutical applications and their requirements.



0.5 m² Stainless Steel Holder

Accomodate an effective filter area of 0.5 m² -2.5 m² up to 5 m²

Ordering Information

	Effective Filtration Area	Filter Plate Material	Inlet/Outlet/Retentate	Seals Materials	Surface Finish	Design Pressure	Industry
H-CFH	<u>01</u>	<u>S</u>	<u>T15</u>	<u>S</u>	<u>A</u>	<u>X</u>	<u>P</u>
	01 0.1 m²	S 316L	T15 DN15	S Silicone	A Mirror Polish	X 0.6Mpa	P Pharmaceutical
	05 0.5 m²		T20 DN20				
	25 2.5 m²		T38 DN38				
	...						

TFF System

Cobetter is able to design TFF system based on customer specifications and requirements especially for Pilot Labs and Process Systems



Pilot System

Throughput (50-100L)

